

RAB12 PAB

产品名称: Rab12 兔多抗

货号: 21076

基因符号: RAB12

描述: Rab12 兔多抗

背景: Rab12 protein is a member of the Rab subfamily of Ras-related small GTPases. Rab12 might regulate vesicular trafficking.

免疫原: A synthetic peptide from the internal region of Rab12

经过测试的应用: ELISA, WB, IHC, IP

推荐稀释比:

ELISA: 1:100-1:5000

WB: 1:200-1:2000

IHC: 1:10-1:100

Concentration: 0.2 mg/ml

种属反应性: Rabbit

形式: Liquid

克隆性: Polyclonal

亚型: IgG


纯化: Purified from serum

成分: PBS (without Mg^{2+} and Ca^{2+}), pH 7.4, 150 mM NaCl, 50% glycerol

种属反应性: Rab12 antibody recognizes Rab12 from vertebrates.

储存条件: Store at $-20^{\circ}C$. Avoid repeated freezing and thawing

蛋白质印迹:



HEK 293T cell lysate was used for Western blot with anti Rab12 兔多抗 (cat # 21076).