

## RAB5-GTP

产品名称: Rab5-GTP 小鼠单抗

货号: 26911

规格: 30  $\mu$ L

基因符号: RAB5

描述: Rab5-GTP 单抗

背景: Rab5 is a member of Rab-family GTPases. Rab5 is a rate-limiting component of the machinery regulating the kinetics of membrane traffic in the early endocytic pathway.

免疫原: Recombinant full length protein of active Gao

经过测试的应用: IP, IHC and IF (**Not applicable for WB since WB denatures the GTPase**)

推荐稀释比: IP: 1  $\mu$ g for 1~2 mg total cellular proteins IHC, IF: 1:50-1:250

**Concentration:** 1 mg/ml

种属反应性: Mouse

形式: Liquid

克隆性: Monoclonal

亚型: IgG1


纯化: Purified from ascites

**Preservative:** No

成分: PBS (without  $Mg^{2+}$  and  $Ca^{2+}$ ), pH 7.4, 150 mM NaCl, 50% glycerol

种属反应性: active Rab5 antibody recognizes active Rab5 proteins of vertebrates.

储存条件: Store at  $-20^{\circ}C$ . Avoid repeated freezing and thawing



免疫沉淀/蛋白质印迹:

### Rab5 Activation Assay.

Purified full-length Rab5A proteins were immunoprecipitated after treated with GDP (lane 1) or GTPγS (lane 2). Immunoprecipitation was done with the anti-Rab5-GTP monoclonal antibody (货号: 26911). Immunoblot was with an anti-Rab5 兔多抗 (货号: 21151).